
ic0943

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***PLEASE DO NOT RETURN  
THIS PRODUCT TO THE STORE***

***FOR QUESTIONS OR ASSISTANCE ON ASSEMBLY, PARTS OR OPERATION,  
CONTACT OUR CUSTOMER SERVICE!***

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***1 866 967.7333 ext. 227***

Thank you for purchasing this product. Please read this instruction manual carefully before operating the hand massager.

Please pay attention to the important safety warnings and keep this manual for future use.

**Note:** EBI reserves the right to revise the design and description of this product without notice.

## CONTENTS

1. Safety warnings.....	2
2. Environment.....	2
3. Restrictions.....	2
4. Maintenance.....	3
5. Troubleshooting.....	3
6. Features and benefits.....	3
7. Components .....	4
8. Operation.....	5
9. Product specifications .....	5
10. Warranty .....	6
FRANÇAIS .....	7

## 1. SAFETY WARNINGS

- This product is not recommended for children.
- Disconnect the massage unit before cleaning.
- Please operate this product according to the instructions in this manual.
- This product is designed for indoor use only.
- Do not exert excessive pressure on the massager. This may significantly slow down or stop the motor, potentially cause overheating and subsequent damage to the massager.
- Press the "stop" button to terminate the massager if discomfort is experienced.
- The total recommended usage time of the massager is 10 minutes.
- Please avoid use if the massage unit is damaged and call your iComfort technician or dealer.
- Do not use wet hands to pull the plug.
- Do not spill liquids onto this product to avoid electric shock or cause damage to this product.
- In the event of a power failure disconnect the massager.

## 2. ENVIRONMENT

The iComfort massager should be used only in well ventilated, normally dry and dust free environments.

## 3. RESTRICTIONS

People with osteoporosis, heart disease and heart pacemakers, pregnancies, injuries or any other health related issues should consult their doctor before using this massage unit.

**Stop immediately the use of the massager if you feel unwell and consult your doctor.**

#### 4. MAINTENANCE

- This product should only be maintained by iComfort certified technicians; otherwise, the warranty will be void.
- Make sure the wall outlet is in good working order before use.
- If the massager will not be used for an extended period of time, store in a dry and dustless environment.
- Avoid product exposure to high temperatures and direct sunlight.
- Clean this product with a dry cloth. Do not use thinner, benzene or alcohol.
- We recommend to use the massager intermittently and not for long time.

#### 5. TROUBLESHOOTING

- It is normal for the motor to make noise while in use.
- If the massager doesn't work normally, please check that the plug and socket are properly connected and the switch is turned on.
- After several use, the massager might power-off because of overheat. It can be used again after half an hour of cooling time.

#### 6. FEATURES AND BENEFITS

Good things really do come in small packages. Fabulously fun and pleasantly powerful, the IC0943 iComfort® hand held massager uses high frequency flap and vibration massage actions. The tapping and vibration movements are done by the 2 rotary heads. With the unique design of the accessory massage heads, you can choose your favorite massage and enjoy its benefits. Great after workouts, long days on your feet, gardening and more. The ergonomic design with long handle makes it easy to target your sore spots even between the shoulder blades. You can appreciate this amazing hand held massager anywhere you desire. Simply set the speed/intensity you want and place it wherever you need relief.

- 2 interchangeable and one permanent massage head system
- With the unique design of the three types of massage heads, you can choose your favorite massage and enjoy its benefits.

## 7. COMPONENTS

1. Speed regulating switch
2. Push switch (On and Off)
3. Massage head
4. Anti-slip handle
5. Massage head accessories

## 8. OPERATION

- Connect the massager to an AC wall outlet.
- Push the 'On/Off' button on the unit's control panel to turn the unit on.
- Set the kneading 'intensity' level by moving the intensity button up or down.
- Press on the 'On/Off' button to turn the unit off.
- Disconnect the massager after use.

The massage strength can be adjusted based on the contact pressure with your body. Please unplug it for cooling if it has been in operation for more than 15 minutes before continuing.

The ic0943 can be used on different body parts.

For the neck and shoulder

For the legs

For the abdomen

You can also enjoy its relaxing effects on other body parts!

## 9. PRODUCT SPECIFICATIONS

Input Voltage: AC110 V - 60 Hz

Power consumption: 25 Watt

## 10. WARRANTY

### LIMITED ONE (1) YEAR WARRANTY

This product is guaranteed against all manufactured faults and defects for a period of one (1) year from the date of purchase.

For repairs or product defects, please call our customer service department at:  
**1-866-967-7333.**

Please have your proof of purchase on hand. If we are not able to replace your defective part, we will send you a service technician for repairs.

**Important notes:** Any damages caused by an accident or an abusive usage are NOT covered under warranty. We also decline all responsibility concerning indirect or consecutive damages. In certain Canadian provinces, it is not permitted to limit indirect or consecutive damages. Therefore it is possible these limitations are not applicable to you. This warranty gives you specific legal rights. You may also have other legal rights which vary from one province to another.

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**Warranty must be active  
within 30 days of purchase.**

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Pictures are for reference only, colors and models may vary without notice. Made in China. iComfort™ is a trademark of EBI Inc. Photographies pour référence seulement. Les couleurs et modèles sont sujet à varier sans préavis. Fabriqué en Chine iComfort™ est une marque de commerce de EBI inc